

KIT FOR THE DETECTION OF METHYLATION OF PROMOTER OF GENES INVOLVED IN COLORECTAL CANCER

AMPLI-SET Colon Cancer Cat. n.1.413

The methylation of the residues of cytosine in the “CpG islands” is very important for the regulation of the genic expression. The hyper- methylation of the “CpG islands” in the promoter region of a gene suppresses the transcription of the same gene. In many tumors the hyper- methylation of the promoter of the suppressor genes, as p16, p15, E-cadherine and other genes as “DAP-kinase”, inhibitor gene of the metastatic progression , 06-metilguanina DNA methyltransferase (MGMT), gene involved in the repair of DNA, Glutathione-S-transferasi (GSPT1) etc.

Plasma and serum of patients carrier of malignant neoplasia contains much genomic DNA than the control subjects. The principle of the assay is the extraction of genomic DNA from plasma or serum , the treatment with bisulfite sodium in order to transform the unmethylated residue of cytosine in uracil, the PCR amplification with specific oligonucleotides for the methylated sequences and unmethylated (MSP:methylation specific PCR) followed by detection on agarose gel. The assessment of the state of hypermethylation of a gene is an appreciable marker of the risk, and allows a precocious diagnosis and a prognosis of a neoplastic diseases.

The kit allows the detection of the methylation of the promoter of the tumour suppressor gene p16 , of the 06-methylguanina DNA methyltransferase gene (MGMT), and of the hMLH1 gene, involved in the repair of “DNA mismatch”, which show respectively a hypermethylation state of the 25-35%, of 40% and of 20-30% in sporadic colorectal cancer.

Principle of Assay: A) extraction of genomic DNA B) Modification treatment with sodium bisulfite C) amplification with specific primers for methylated and unmethylated sequences of the promoter of the genes: p16 - hLMH1 – MGMT, D) detection on agarose gel.

Applicability: On extracted and purified genomic DNA from whole blood samples.

Numbers of Tests: 24

MODIFICATION	
NaOH	-20°C
Reagent A (sodium bisulfite)	-20°C
Reagent B (idrochinone)	-20°C
Reagent C (resin)	-20°C
Reagente (carrier)	-20°C
Diluyente	-20°C
Sterile H ₂ O	-20°C
AMPLIFICATION	
Mix PCR METHYLATED MGMT	-20°C
Mix PCR UNMETHYLATED MGMT	-20°C
Mix PCR METHYLATED p16	-20°C
Mix PCR UNMETHYLATED p16	-20°C
Mix PCR METHYLATED hLMH1	-20°C
Mix PCR UNMETHYLATED hLMH1	-20°C
sterile H ₂ O	-20°C
Taq Polymerase (5U/□l)	-20°C
Unmethylated DNA Control	-20°C
Methylated DNA Control	-20°C

REAGENTS and STORAGE

Stability: over 12 months if correctly stored.

References:

- 1) sanchez-cespedes m. Et al cancer res 60, 892-895 (2000)
- 2) herman j. G. Et al proc. Natl. Acad. Sci. Usa 95, 6870-6875 (1998)
- 3) esteller m. Et al oncogene 16, 2413-2417 (1998)
- 4) esteller m. Et al cancer res 59, 67-70 (1999)
- 5) esteller m. Et al cancer res 59, 793-797 (1999)
- 6) leon s. A. Et al. Cancer res 37, 646-650 (1977)
- 7) stroun m. Et al oncology 46, 318-322 (1989)
- 8) shapiro b. Et al cancer 51, 2116-2120 (1983)
- 9) wong i. H. N. Et al cancer res 59, 71-73 (1999)
- 10) baylin s. B. Adv. Cancer res. 72, 141-196 (1998)
- 11) belinsky s. A. Proc. Natl. Acad. Sci. Usa 95, 11891-11896 (1998)


ANALYSIS OF RESULTS

Every sample of DNA, after treatment with sodium bisulfite, has been amplified both with METHYLATED (M) mix and UNMETHYLATED (UM) mix using specific oligonucleotides included in the kit for the methylated and unmethylated sequences of the p16, hLMH1 e MGMT genes

Every PCR product must be run on agarose gel separately (it is suggest to put the product of amplification of every sample near each other on the gel). Samples where a methylated state of the promotor of the gene is present , give a band both with the METHYLATED mix and the UNMETHYLATED mix.

Samples where an methylated state of the promotor of the gene isn't present, give a band only with the UNMETHYLATED mix. Specifically, the analysis of the status of methylation of the promotor of the p16 gene gives a band of 150 bp with the METHYLATED MIX and a band of 151 bp with UNMETHYLATED mix . The analysis of the status of methylation of the promotor of MGMT gene gives a band of 98 bp with METHYLATED mix and of 106 bp with the UNMETHYLATED mix PCR. The analysis of the status of methylation of the promotor of the hMLH1 gene gives a band of 105 bp with the METHYLATED mix PCR and a band of 114 bp with the UNMETHYLATED mix PCR

[Example: analysis of methylation status of p16 gene](#)



N.B. The DNA unmethylated Control , after treatment with sodium bisulfite, will show a pCR product only with the **Unmethylated Mix.**