





KIT FOR THE DETECTION OF 4G/5G POLYMORPHISM IN THE PROMOTER OF THE GENE PAI-1

AMPLI-SET-PAI-4G/5G

Cat. n.1.351

The development of thrombotic pathologies is one of the major cause of morbidity and mortality. The alteration of haemostatic system is the central mechanism of the thrombotic events. The causes of this alterations may be also genetic. It underlines the importance of the interaction between genes and environment in the thrombotic pathologies. The inhibitor of the activator of Plasminogen (PAI-1) is the main inhibitor of the activators of plasminogen (t-PA and u-PA). An increase or a decrease of PAI-1 can be responsible of thrombotic or haemorrhagic events respectively. Recently, the detection of a polymorphism in the promoter region of the PAI-1 gene has lead to evaluate PAI-1 as a factor risk of venous thrombosis. Particularly, the insertion/deletion of a residue of Guanosine in the region 675 of the promoter of the gene, named 4G/5G, has been related with an increase of the activity of PAI-1. This polymorphism involves the bond of nuclear proteins responsible of the transcription of the gene PAI-1. The site 4G may links only an "enhancer" region, whereas the allele 5G links both the "enhancer" region and the "suppressor" one, resulting in a lower level of transcription. "In vivo" studies show that levels of PAI-1 are 25% higher in homozygotic subject for the allele 4G than in homozygote subject for allele 5G. Homozygosis for the allele 5G is present in 25% of normal population, and in 10% of pathological population for thrombotic events. The relationship between polymorphism 4G and the risk of venous thrombosis needs of further investigation.

Probably, the presence of the homozygosis 4G becomes a risk factor when is joined to other genetic or environmental risk factors.

The detection of 4G/5G PAI-1 polymorphism is performed with the amplification with specific primers of a fragment of 228 bp, followed by restriction section due to *Eci I* enzyme. The 5G polymorphism introduces a restriction site for the enzyme. The product of PCR of the allele 4G isn't cut, whereas the allele 5G produces two fragments of 188bp and 40 bp.

Principle of method: A) extraction of genomic DNA B) amplification C) enzymatic digestion D)detection on agarose gel.

Applicability: On extracted and purified genomic DNA from whole blood samples.

Tests: 45

REAGENTS AND STORAGE


AMPLIFICATION	
PCR mix PAI 4G/5G	-20°C
H ₂ O sterile	-20°C
Taq Polymerase (5U/ l)	-20°C
Eci I enzyme (2U// 1)	-20°C
Digestion buffer 10X	-20°C
BSA 100X	-20°C
Positive Control 4G homozigous	-20°C
Positive control 5G homozigous	-20°C

Stability: over 12 months if correctly stored.

References:

- 1. Hemostasis and Thrombosis: Basic Principles and Clinical Practice, 4th ed. Philadelphia, Pa: 2001: 975-1002.
- 2. Arterioscler Thromb. 1991; 11: 183-190.
- 3. Thromb Haemost. 1995; 74: 837-841.
- 4. Proc Natl Acad Sci USA 1995; 92: 1851-1855.
- 5. T hromb Haemost. 1998; 79:354-358.
- 6. Arterioscler Thromb Vasc Biol.
- 7. Arch Pathol Lab Med 2002; 126: 1401-1404.

ANALYSIS OF RESULTS



 $\mathbf{M} = \text{Marker 50bp ladder}$

- 1) = Heterozygote 4G/5G
- 2) = Homozygote 4G (control DNA)
- 3) = Homozygote 5G (control DNA)

Heterozygote	Homozygote	Homozygote
4G/5G	4G/4G	5G/5G
228 bp	228 bp	188 bp
188 bp	•	40 bp
40 bp		

Usually, the resolution of agarose gel weakens the visualization of the band of 40 bp. The molecular diagnosis is guaranteed from the visualization of other fragments.