

HPV SCREENING

Amplify-HPV Screening

Cat. n. 1.111

The infection by HPV (Human Papilloma Virus) is one of the main cause of cervical intraepithelial neoplasia and of the ano-genital carcinoma. Among the different types of HPV that infect the genital and anal mucosa, the 16, 18 and 33 have been detected mainly in invasive cervical carcinoma ; the 6 and 11 are been related mainly to benign neoplasia.

The Amplify-Set HPV screening allows to detect, using the PCR, Polymerase Chain Reaction, the viral DNA of HPV of these types: 6, 11, 13, 16, 18, 26, 31, 33, 35, 39, 42, 43, 44, 45, 51, 52, 53, 54, 56, 57, 58, 59 in biological samples as omogenate tissue and cells preparations. The amplified region by the primers used in the kit is the L1 region of the viral genome. Analytical sensitivity is greater than 16 viral genomes / μ l, specificity 100 %.

Principles of assay: A) Extraction of genomic DNA B) amplification C) detection on agarose gel.

Applicability: on extracted and purified DNA of omogenate tissue and cells preparations.

Tests: 40.

REAGENTS AND STORAGE


AMPLIFICATION

PCR mix *	-20°C
H ₂ O sterile	-20°C
Taq Polymerase (5U/ μ l)	-20°C
Positive control	-20°C

* Mix PCR contains an internal PCR control that produces a fragment of 800 bp.

Stability: over 12 months if correctly stored.

ANALYSIS OF RESULTS



The amplification yield of DNA-HPV positive samples is of 420 bp (lane 1). The amplification yield of negative samples is of 800 bp (internal control of amplification) (lane 2)

References:

- Journal of Clinical Microbiology* 33, 690-695 (1995).
- Journal of Virology* 69, 3074-3083 (1995)
- Journal of Investigative Dermatology* 105, 367-371 (1995)
- Virology* 184, 492-503 (1991)