

## KIT FOR THE DETECTION OF C677T POLYMORPHISM OF THE METHYLENETETRAHYDROFOLATE REDUCTASE (MTHFR) GENE

### AMPLI-SET-MTHFR C677T Cat. n. 1.300

The increase of the level of homocysteine may depend on a metabolic block of the transformation of homocysteine in cystathionine or on the unsuccessful methylation of homocysteine in methionine. The enzyme MTHFR catalyses the reduction of 5,10-methylenetetrahydrofolate in 5-methyltetrahydrofolate, the predominant form of circulating folate and donor of carbon in the process of re-methylation of homocysteine in methionine.

A mutation C-T, which inserts a Valine instead of a Alanine, is associated with a reduced activity and an increased thermolability of this enzyme. Homozygote subjects for the mutation show a significative increase of the plasmatic level of circulating homocysteine, due to the unsuccessful conversion in methionine.

Fasting hyper-homocysteinemia (increased plasmatic level of circulating homocysteine) is associated to an increased risk of vascular cerebral, peripheral and coronary diseases

The detection of the MTHFR (C-T) is carried out starting with an amplification using specific primers of a fragment 198 bp, following by a restriction section due to Hinf I enzyme. The mutation is confirmed by the detection of a restriction cleavage for the Hinf I enzyme. So, the amplification product of the normal allele isn't cut, whereas the one of the mutant allele produces two fragments of 175bp and 23bp.

**Principle of Assay:** A) extraction of genomic DNA B) amplification C) enzymatic digestion D) detection on agarose gel.

**Applicability:** On extracted and purified genomic DNA from whole blood samples.

**Tests:** 45.

### ANALYSIS OF RESULTS

The yield of amplification is a fragment of 198 bp. the PCR fragment containing the mutation is cleaved into two fragments (175 and 23 bp).

### REAGENTS AND STORAGE

EXTRACTION	
Reagent 1 (Separation buffer)	+4°C
Reagent 2 (Lysis buffer)	+4°C
Reagent 3 (Deproteination solution)	+4°C
Reagent 4 (Proteinase k)	-20°C
Reagent 5 (NaCl solution)	+4°C
AMPLIFICATION and DIGESTION	
PCR mix	-20°C
H <sub>2</sub> O sterile	-20°C
Taq Polymerase (5U/μl)	-20°C
Hinf I enzyme (10U/μl)	-20°C
Digestion buffer 10X	-20°C
Positive control (Homozygote mutated)	-20°C
DETECTION	
Agarose gel 4% for electrophoresis in TBE 1X	R.T.
Loading buffer 10 X	R.T.
TBE 10 X	R.T.
Marker ladder 100 bp	-20°C


**Stability:** over 12 months if correctly stored (Agarose gels, if protected by light, can be stored 1 year at room temperature).

**Required Equipment:** tubes 1.5 ml screw caps; tubes 15 ml in polypropylene refreshed tube-tray; plugged (aerosol-barrier) sterile tips; glass Pasteur; tubes NO-OIL

**Reagents not included with the kit:** bidistilled sterile water, chloroform, Ethanol 70%, Ethanol 100%.

**Required instruments:** microcentrifuge with fixed rotor for 2 ml tubes (13-15.000 g); Pre-PCR and Post-PCR Pipettes (0.5 – 20 μl, 10 – 100 μl, 20 – 200 μl, 200 – 1000 μl); Programmable Thermalcycler; Class II Bio Hazard Biological Safety Cabinet; Electrophoresis system with power supply; UV Transilluminator; Photodocumentation System.

*All materials requested must be DNase/RNase free, sterile and disposable. It is suggested to operate under safety cabinet and in ice. It is absolutely necessary to use gloves during the handling of samples and reagents.*



### LEGENDA GEL

- 1) Amplification product of a DNA HOMOZYGOUS RECESSIVE subject
- 2) Restriction cleavage with HINF I of the sample 1
- 3) Amplification product from DNA ETEROZYGOSIS subject
- 4) Restriction cleavage with con HINF I of the sample 2
- 5) Restriction cleavage on eterozygote control
- 6) restriction cleavage on HOMOZYGOUS NORMAL CONTROL.
- 7) Marker VIII

### REFERENCES:

- Frost P et al. *Nat Genet* 1995 May;10(1):111-3.  
Margaglione M et al. *Thromb Haemost* 1998 May;79(5):907-11.